Reference-based PCR Duplicate Removal

Review

November 27, 2017

What does a PCR-duplicate look like?

Same alignment position

• Chromosome RNAME (SAM col 3)

Position
 Pos (SAM col 4)

• Strand (strand specific?) FLAG (SAM col 2)

• Soft Clipping CIGAR (SAM col 6)

• Same Unique Molecular Index QNAME (SAM col 1) (UMI or "randomer")

• Single-end vs Paired-end? FLAG (SAM col 2)

SAM: alignment section mandatory fields

• SAM format: https://samtools.github.io/hts-specs/SAMv1.pdf

| Col | Field | Туре | Description | |
|-----|--------------|--------|---------------------------------------|--|
| 1 | QNAME | String | Query template NAME | |
| 2 | FLAG | Int | bitwise FLAG | |
| 3 | RNAME | String | Reference sequence NAME | |
| 4 | POS | Int | 1-based leftmost mapping POSition | |
| 5 | MAPQ | Int | MAPping quality | |
| 6 | CIGAR | String | CIGAR string | |
| 7 | RNEXT | String | Ref. name of the mate/next read | |
| 8 | PNEXT | Int | Position of the mate/next read | |
| 9 | TLEN | Int | observed Template LENgth | |
| 10 | SEQ | String | segment SEQuence | |
| 11 | QUAL | String | ASCII of Phred-scaled base QUALity+33 | |

SAM: parsing the bitwise FLAG

https://broadinstitute.github.io/picard/explain-flags.html

```
if ((flag & 4) != 4):
    mapped = True
```

| Bit | | Description | |
|-------------|-------|--|--|
| 1 | 0x1 | template having multiple segments in sequencing | |
| 2 | 0x2 | each segment properly aligned according to the aligner | |
| 4 | 0x4 | segment unmapped | |
| 8 | 0x8 | next segment in the template unmapped | |
| | 0x10 | SEQ being reverse complemented | |
| 32 | 0x20 | SEQ of the next segment in the template being reverse complemented | |
| 64 | 0x40 | the first segment in the template | |
| 128 | 0x80 | the last segment in the template | |
| 256 | 0x100 | secondary alignment | |
| 512 | 0x200 | not passing filters, such as platform/vendor quality controls | |
| 1024 | 0x400 | PCR or optical duplicate | |
| 2048 | 0x800 | supplementary alignment | |



Soft clipping

- What is it?
- What does it look like?
- Why would something be soft clipped?
 - Sequence error/heterozygosity
 - Over-penalizing indels
 - Splicing with just a few nucleotides in an exon
 - Novel splicing
- Where in the alignment could soft clipping occur?

How do you know if your sequence was soft clipped? The CIGAR string!

CIGAR: CIGAR string. The CIGAR operations are given in the following table (set '*' if unavailable):

| Op | BAM | Description |
|----|-----|---|
| M | 0 | alignment match (can be a sequence match or mismatch) |
| I | 1 | insertion to the reference |
| D | 2 | deletion from the reference |
| N | 3 | skipped region from the reference |
| S | 4 | soft clipping (clipped sequences present in SEQ) |
| H | 5 | hard clipping (clipped sequences NOT present in SEQ) |
| P | 6 | padding (silent deletion from padded reference) |
| = | 7 | sequence match |
| X | 8 | sequence mismatch |

```
10M 5M1024N5M

...ATTGTCCATT... ...ATTGTAGT...GCCCATT...
10M ATTGACCATT 2S8M ATTGT CCATT

...ATTGTCCATT... ...ATTGTCCATT...
GGTGTCCATT...
GGTGTCCATT
```

Example CIGAR strings from real data

| 1527M113N38M | 13M1185N53M | 23M686N43M | 40M1I25M |
|--------------|------------------|----------------|--------------|
| 1536M284N29M | 13M83N53M | 25M246N41M | 43M32325N23M |
| 1S45M301N20M | 13M85N53M | 25M470N41M | 44M5064N22M |
| 1561M4S | 14M903N52M | 27M2645N39M | 46M2001N20M |
| 1564M1S | 15M470N51M | 28M138N38M | 49M211N17M |
| 2S31M745N33M | 16M1540N50M | 28M1794N38M | 51M97N15M |
| 2562M2S | 17M1097N49M | 28M3349N38M | 52M1198N14M |
| 2563M1S | 18M12872N48M | 29M2D37M | 52M1904N14M |
| 2564M | 19M20545N47M | 30M1244N36M | 54M12S |
| 3S18M104N45M | 20M2979N46M | 31M1043N35M | 54M407N12M |
| 3S63M | 20M456N46M | 31M271N35M | 55M954N11M |
| 4562M | 20M631N46M | 35M113N31M | 56M10S |
| 5S61M | 21M982N45M | 35M128N31M | 57M5055N9M |
| 6542M284N18M | 22M103N44M | 35M289N31M | 58M8S |
| 6560M | 22M138N44M | 35M5284N31M | 62M4S |
| 7S59M | 22M4529N44M | 36M103N30M | 63M3S |
| 11M4797N55M | 23M1290N25M2D18M | 36M12071N29M1S | 65M1S |
| 11M592N55M | 23M15018N43M | 36M1496N30M | 66M |

Your algorithm!

Due January 12, 2018

Given a SAM file of uniquely mapped reads, remove all PCR duplicates (retain only a single copy of each read)

- Samtools sort
- Adjust for soft clipping
- Single-end reads
- Known UMIs
- Considerations:
 - Millions of reads avoid loading everything into memory!
 - Be sure to utilize functions appropriately
 - Appropriately commend code and include doc strings
- CHALLENGE: Include options for
 - Single-end vs paired-end
 - Known UMIs vs randomers
 - Choice of duplicate written to file

You MUST

- Write Python 3 compatible code
- Include the following argparse options

```
    -f, --file required arg absolute file path
    -p, --paired optional arg designates file is paired end (not single-end)
    -u, --umi optional arg designates file containing the list of UMIs (unset if randomers instead of UMIs)
    -h, --help optional arg docs)
```

- If your script is not capable of dealing with a particular option (ex: no paired-end functionality), your script should print an error message and quit
- Output the first read encountered if duplicates are found
 - You may include an additional argument to designate output of a different read (highest quality or random or ???)
- Output a properly formatted SAM file with "_deduped" appended to the filename
- Name your python script < your_last_name > _deduper.py